

DEMOPHORIOUS IS DEDICATED TO MANUFACTURING MEDICAL RELATED PRODUCTS, OF THE HIGHEST QUALITY AND STANDARDS. DEMOTEK **BLOOD BAGS** ARE MANUFACTURED UNDER STRICT QUALITY CONTROL TO ENSURE SAFETY, RELIABILITY AS WELL AS EXCELLENT PRODUCT QUALITY. D-VAC BLOOD COLLECTION SYSTEM HAS BEEN DESIGNED TO CONFORM WITH INTERNATIONAL STANDARDS AND SAFETY. D-TEK SUTURES HAS ACHIEVED TO BE BRANDED AMONG THE TOP BRANDS DUE TO THE COMPANY'S COMMITMENT OF MANUFACTURING HIGH QUALITY SUTURES. DEMOTEK MEDICAL DISPOSABLES HAVE PROVED TO BE ONE OF THE MOST REPUTABLE BRANDS IN GENERAL HOSPITAL SUPPLIES.

DEMOTEK BLOOD BAGS ARE MANUFACTURED UNDER STRICT QUALITY CONTROL TO ENSURE SAFETY, RELIABILITY AS WELL AS EXCELLENT PRODUCT QUALITY.

**COMMITTED TO PRODUCT DEVELOPMENT
QUALITY & RELIABILITY.**

**DEMOTEK BLOOD BAGS ARE CONTINUOUSLY
UPGRADED TO MEET THE NEW CHALLENGING
DEMANDS IN MEDICAL CARE** AS WELL AS THE
GROWING NEEDS FOR SAFETY ENHANCEMENT
IN BLOOD RELATED PRODUCTS LIKE THE
PREVENTION OF ADVERSE EFFECTS ON BLOOD
TRANSFUSION.

DEMOTEK HIGH QUALITY BLOOD BAGS ARE
SUITABLE FOR THE COLLECTION OF BLOOD AND
BLOOD COMPONENTS AND FOR THE FILTRATION
OF LEUKOCYTES FROM WHOLE BLOOD, RED
CELLS AND PLATELET PRODUCTS.

BLOOD BAGS

ARE MANUFACTURED UNDER STRICT
QUALITY CONTROL TO ENSURE
SAFETY RELIABILITY & EXCELLENT
PRODUCT QUALITY.

BB

TRIPLE BLOOD BAGS

SUITABLE FOR PRESERVATION AND TRANSFER OF HUMAN BLOOD OR ITS COMPONENTS.

COLLECTION & SEPARATION OF WHOLE BLOOD:

3X 	BLOOD COMPONENTS: <ul style="list-style-type: none"> • RED BLOOD CELLS • PLATELETS • PLASMA 	OBTAINED THROUGH THE PROCESS OF: <ul style="list-style-type: none"> • CENTRIFUGATION • EXTRACTION 	350ml
			450ml

A PLATELET STORAGE BAG CAN BE INCORPORATED FOR **5-DAYS** PLATELET STORAGE.



TRIPLE BAG COMES WITH: • CPDA-1 • CPD/SAG-M SOLUTION

LABEL REGULATORY

ISO3826-1
A

ISO15223-1
B

ISBT128
C


LR

demotek

NEEDLE: 16G

EACH 100ML OF CPDA-1 CONTAINS:

- Citric Acid (Anhydrous) 0.299g.
- Monobasic Sodium Phosphate (Monohydrate) 0.222g.
- Adenine 0.0275g.
- Sodium Citrate (Dihydrate) 2.63g.
- Dextrose (Monohydrate) 3.19g.
- Water for injection to 100ml.



CAUTIONS:

1. Do not use the bag more than 15 days after opening the aluminum foil pack.
2. Do not use if there is any visible sign of deterioration.

5 UNITS

STORAGE:

1. Unused pack may be stored at room temperature less than 25oC, relative humidity less than 90% and pack containing blood should be kept between +2oC and 6oC.
2. Avoid excessive heat and direct sunlight.



REF BBS090



1D30BBS090



Demophorus Limited
196 Arch. Makariou III Ave.,
3030 Limassol, Cyprus, EU.
info@demophorus.com

LOT 111111




0000111111

LABEL REGULATORY

A ISO3826-1
PLASTICS COLLAPSIBLE CONTAINERS FOR HUMAN BLOOD AND BLOOD COMPONENTS.
PART 1: CONVENTIONAL CONTAINERS.

B ISO15223-1
MEDICAL DEVICES. SYMBOLS TO BE USED WITH MEDICAL DEVICE LABELS, LABELLING AND INFORMATION TO BE SUPPLIED. **PART 1: GENERAL REQUIREMENTS.**

C ISBT128
STANDARD BARCODING SYSTEM.



demotek TRIPLE BLOOD BAG 450ML

For blood collection volume 450ml
CPD 63ml

- 600ml
- CPD 63ml
- 450ml

Each 100ml CPD contains:
Citrates 2.5g
Dextrose 0.5g
Mannitol 0.5g
Sodium Citrate 1.5g
Sodium Citrate Dihydrate 1.5g
Sodium Citrate Trihydrate 1.5g
Sodium Citrate Tetrahydrate 1.5g
Sodium Citrate Pentahydrate 1.5g
Water for Injection 100ml

The bag system is suitable for preservation or transfer of human blood or its components

Needle: 18G
Store unused bag at 15°C-25°C
Store blood at 4°C±2°C
Transfer set must have a filter
Rx ONLY

CE 0123
2013.10.17
STERILE
DEHP

Donation No.
Collection Date Expiry Date
ABO/Rhd Group
Centre Identity

Demotek Health Care
199 Athinikou 10 Ave.
5700 Limassol Cyprus EU
info@demotek.com

2016.10.16

LOT 131017
00001 31017

REF BBTB4
1DE08 BTTB4

Do Not Fix ONLY
CE 0123
2013.10.17
BBTTB4
31017



demophorius

HEALTHCARE

DEMOPHORUS.COM

196 ARCHBISHOP MAKARIOS III AVE
LIMASSOL 3030, CYPRUS

E. INFO@DEMOPHORUS.COM

T. +357 25 74 98 96

F. +357 25 74 99 07